

RADIESSE® for Hands

Smooth Over Hand Wrinkles

Just like our faces, our hands change over time. Veins and tendons become more prominent, and the youthful appearance fades.

RADIESSE® is the first dermal filler proven to immediately restore volume loss in the back of hands¹, giving you smooth, natural-looking, and long-lasting results.²



HOW RADIESSE CAN ENHANCE YOUR HANDS

- [Restoring Lost Volume](#)
- [Minimizing the Visibility of Veins and Tendons](#)
- [Providing a Smooth, Natural-Looking, and Long-Lasting Result²](#)
 - [Opaque Gel Carrier²](#)



- **Restoring Lost Volume**
Your hands can show the signs of aging through volume loss.

Immediate results **that last**

The hand and facial filler works to add volume under the skin, and over time, the benefits of RADIESSE continue by stimulating your body's own natural collagen.

Zondra, 49

RESULTS AFTER 2 WEEKS*



*Treatment volume: 1 syringe (1.5 cc) per hand. Individual results may vary. Unretouched photos of actual patient.

HANDS BEFORE & AFTER

What are RADIESSE® and RADIESSE® (+)?

RADIESSE® and RADIESSE® (+) are dermal fillers that are used for smoothing moderate to severe facial wrinkles and folds, such as nasolabial folds (the creases that extend from the corner of your nose to the corner of your mouth). RADIESSE® is also used for correcting volume loss in the back of the hands.

RADIESSE® and RADIESSE® (+) IMPORTANT SAFETY INFORMATION

Who should not use RADIESSE® or RADIESSE® (+)?

You should not use RADIESSE® or RADIESSE® (+) if you have an allergy to any component of the product, if you have a history of severe allergies, if you have a bleeding disorder, or if you are pregnant or breastfeeding. You should not use RADIESSE® (+) if you have an allergy to lidocaine or medicines like it.

What is the most important information I should know about RADIESSE® and RADIESSE® (+)?

One of the risks with using these products is unintentional injection into a blood vessel. The chances of this happening are very small, but if it does happen, the complications can be serious, and may be permanent. These complications, which have been reported for facial injections, can include vision abnormalities, blindness, stroke, temporary scabs, or permanent scarring of the skin. If you have changes in your vision, signs of a stroke (including sudden difficulty speaking, numbness or weakness in your face, arms, or legs, difficulty walking, face drooping, severe headache, dizziness, or confusion), white appearance of the skin, or unusual pain during or shortly after treatment, you should notify your health care practitioner immediately.

As with all procedures that involve an injection through the skin, there is a risk of infection.

Do not use RADIESSE® or RADIESSE® (+) if you have a skin infection until it has healed.

It is not known if RADIESSE® or RADIESSE® (+) is safe or effective in the lips, or in the area around the eyes.

The microspheres in RADIESSE® and RADIESSE® (+) can be seen in X-rays and CT Scans. It is very important that you tell your health care provider that you have had RADIESSE® or RADIESSE® (+) dermal filler.

If you have a history of herpes, you may experience a herpes breakout after receiving RADIESSE® or RADIESSE® (+).

Injection in the back of the hand may result in temporary difficulty performing activities. RADIESSE® may cause nodules, bumps or lumps in the back of the hand and can last up to 1 year.

You should minimize strenuous activity and avoid extensive sun or heat exposure for about 24 hours after treatment and until any swelling or redness has resolved.

What should I tell my doctor before using RADIESSE® or RADIESSE® (+)?

Tell your health care provider if you are taking blood thinners or medicines that can interfere with the clotting of blood, such as aspirin or warfarin. These medicines might make it more likely that you will experience bruising or bleeding at the injection site. Tell your health care provider if you have any diseases, injuries or disabilities of the hand, if you have a history forming large, raised scars or if you have had any other skin treatments such as skin peels.

What are the most common adverse events with RADIESSE® or RADIESSE® (+)?

The most common adverse events seen in clinical studies of RADIESSE® used in the hands include bruising, redness, swelling, pain, itching, nodules or bumps/lumps, difficulty performing activities, loss of sensation and other local side effects. The most common adverse events seen in clinical studies of RADIESSE® or RADIESSE® (+) used in the face include bruising, redness, swelling, pain, itching and other local side effects.

These are not all of the possible side effects with RADIESSE® or RADIESSE® (+). Merz collects information about adverse events seen with RADIESSE® and RADIESSE® (+) outside of clinical studies. These events are included in the RADIESSE® and RADIESSE® (+) Patient Information Guide based on an assessment of seriousness and potential causal relationship to RADIESSE® or RADIESSE® (+). Please see the Patient Information Guide available at www.radiesse.com for list of these events. Tell your health care provider about any side effects that bother you or do not go away.

Important: For full safety information, please visit www.Radiesse.com or call MyMerz Solutions at 1-844-469-6379

RADIESSE® and RADIESSE® (+) are available by prescription only.

References:

1. Goldman MP, Moradi A, Gold MH, et al. Calcium hydroxylapatite dermal filler for treatment of dorsal hand volume loss: results from a 12-month, multicenter, randomized, blinded trial. *Dermatol Surg.* 2018;44(1):75-83.
2. RADIESSE [Instructions for Use for the dorsum of the hand], Franksville, WI: Merz North America, Inc.; 2016.
3. Edelson KL. Hand recontouring with calcium hydroxylapatite (Radiesse). *J Cosmet Dermatol.* 2009;8:44-51.
4. Berlin AL, et al. *Dermatol Surg.* 2008;34(suppl 1):S64-S67 and/or Yutskovskaya Y, et al. *J Drugs Dermatol.* 2014;13(9):1047-1052.